

**EUROPE**

LIEGE SCIENCE PARK • 4102 Seraing • Belgium • Tel.: +32 4 372 74 00  
 Fax: +32 4 372 75 00 • info@eurogentec.com • www.eurogentec.com

**NORTH AMERICA**

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA  
 Tel.: 510-791-9560 • Toll-free: 800-452-5530 • Fax: 510-791-9572  
 info.anaspec@eurogentec.com • www.eurogentec.com

# HGS Diamond *Taq*<sup>®</sup> DNA Polymerase

## Specification Sheet

### Reference: TAQ-I011

Eurogentec products are sold for research or laboratory use only and are not to be administered to humans or used for medical diagnostics. For medical diagnostics, please use the TAQ-I010 references.

**Source**

HGS Diamond *Taq*<sup>®</sup> is a highly thermostable enzyme produced and purified from recombinant *Escherichia coli* bacterium containing the *Thermus aquaticus* DNA Polymerase gene.

**Intended use**

HGS Diamond *Taq*<sup>®</sup> lacks *Taq* I restriction endonuclease activity. The enzyme shows very good fidelity and catalyzes 5'→3' polymerization-dependent exonuclease replacement activity with no detectable 3'→5' exonuclease activity. HGS Diamond *Taq*<sup>®</sup> is a chemically modified Hot Start *Taq* DNA polymerase, which completely lacks any activity before activation to avoid non-specific priming at low temperature. This enzyme requires a 10 minutes activation step at 95°C to reach maximal initial activity. During the PCR the rest of its activity is released. It is heat-degraded at a much lower rate as commonly used *Taq* DNA polymerase. DNA fragments as long as 2 kb can be efficiently amplified. HGS Diamond *Taq*<sup>®</sup> DNA polymerase provides efficient amplification of specific products without amplifying non-specific products or primer dimers. HGS Diamond *Taq*<sup>®</sup> is particularly suited for diagnostic PCR & qPCR applications that require high sensitivity and ultra low levels of bacterial & fungal and/or highly specific amplification. The GMP manufacturing & purification processes minimize the risk of false positive results due to residual DNA contamination (bacterial or fungal). The enzyme is QC-tested to verify that < 1fg of genomic *E. coli* DNA (or 0.2 copy) is present in a standard aliquot containing 1 unit of *Taq*. Bioburden is guaranteed ≤ 10 CFU/ml, but is typically = 0 CFU/ml.

**Package contents**

Reference	Units	Volume	Concentration	Volume HGS Diamond <i>Taq</i> <sup>®</sup> reaction buffer (10 X)*	Volume 25 mM MgCl <sub>2</sub>
TAQ-I011-100 (sample)	100	20 µl	5 U/µl	1 ml	1 ml
TAQ-I011-1000	1000	200 µl	5 U/µl	6 ml	6 ml
TAQ-I011-5000	5000	1 ml	5 U/µl	30 ml	30 ml
TAQ-I011-25000	5 x 5000	5 x 1 ml	5 U/µl	5 x 30 ml	5 x 30 ml

\* 150 mM Tris-HCl pH 8.5 (at 19°C), 500 mM KCl and stabilizers.

**Shipping conditions**

Shipping at room temperature

**Storage conditions**

Storage at -20°C is recommended

**Storage and dilution buffer**

20 mM Tris-HCl, 1 mM DTT, 0.1 mM EDTA, 0.1 M KCl, 0.5% (v/v) Nonidet P40, 0.5% (v/v) Tween 20, 50% (v/v) glycerol, pH 8.0 (19°C) and stabilizer.

**Enzyme Specifications**

Each lot of enzyme, buffer and MgCl<sub>2</sub> is functionally tested and quality controlled to ensure the following specifications of the IVD-GMP products.

<b>Appearance</b>	Colourless solution
<b>Identity</b>	MW approx. 95 kDa
<b>Volume activity</b>	≥ 5 Units/µl
<b>Specific activity</b>	NA
<b>Purity</b>	> 98%
<b>Performance test on λ DNA (PCR)</b>	0.5 kb fragment positive down to 5 pg
<b>Performance test on 18S DNA (PCR)</b>	0.1 kb fragment positive down to 10 pg
<b>Absence of Ribonucleases</b>	Not detectable up to 10 U, 1 h, 37°C
<b>Absence of Endonucleases</b>	Not detectable up to 30 U, 16h, 65°C
<b>Absence of Exonucleases</b>	Not detectable up to 30 U, 16h, 65°C

## EUROPE

LIEGE SCIENCE PARK • 4102 Seraing • Belgium • Tel.: +32 4 372 74 00  
Fax: +32 4 372 75 00 • info@eurogentec.com • www.eurogentec.com

## NORTH AMERICA

ANASPEC - 34801 Campus Drive • Fremont, CA 94555 • USA  
Tel.: 510-791-9560 • Toll-free: 800-452-5530 • Fax: 510-791-9572  
info.anaspec@eurogentec.com • www.eurogentec.com



# Eurogentec

Experience true partnership

<b>Absence of Nicking activity</b>	Not detectable up to 30 U, 16h, 65°C
<b>5'-3' exonuclease activity</b>	Positive
<b>Hotstart (SYBR® Green qPCR)</b>	No detectable amplification without a heat activation step. Detectable amplification with an activation step for 10 min at 95°C.
<b>Residual <i>E.coli</i> DNA</b>	< 1fg / Taq Unit
<b>Bioburden</b>	≤10 CFU/ml

## 8. Unit definition

One unit is defined as the amount of enzyme that incorporates, after activation step, 10 nmoles of dNTPs into acid insoluble form in 30 minutes at 74 °C.

### Reaction Conditions

#### For a 100 µl Reaction

HGS Diamond Taq <sup>®</sup> Reaction Buffer (10x)	10 µl
MgCl <sub>2</sub> solution	6 µl (1.5 mM)
HGS Diamond Taq <sup>®</sup>	0.8 to 2.5 units
dNTP	200 µM each dNTP
Primers	0.1 nmol each
H <sub>2</sub> O	As required
DNA template	As required

### Magnesium

*This DNA polymerase is a magnesium-dependent enzyme. We recommend increasing the magnesium concentration for long DNA fragments. Excess Mg<sup>2+</sup> stabilizes the DNA double strand and consequently prevents complete denaturation of DNA, which reduces the extension yield. It may also stabilize spurious primer/template annealing, thus decreasing specificity.*

### Recommendation

*Homogenize HGS Diamond Taq<sup>®</sup> solution by flipping the tube 4 to 5 times.*

## Cycling conditions

### Classical PCR protocol used for 500 bp lambda DNA amplification\*

25 cycles	{	95°C	10 min
		(enzyme activation + DNA denaturation)	
	{	94°C	30 sec
		T <sub>m</sub> -2°C	30 sec
		72°C	1 min/kb
	{	72°C	7 min
		4°C	end temperature

*\*Condition will vary from reaction to reaction and may need optimization for maximal performances. Duration and temperature for denaturation and annealing steps depend on the type of cyclers and primers design.*

## Disclaimer

### Roche Enzyme Products Conveying Polymerase Patent Rights

Use of this product is covered by the following US patent 6,127,155. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim, no right to perform any patented method, and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

### Chemically-Modified Hot-Start Polymerase Reagents and Kits

Use of this product is covered by one or more of the following US patents and corresponding patent claims outside the US: 6,127,155, 5,677,152 (claims 1 to 23 only), 5,773,258 (claims 1 and 6 only), and claims outside the US corresponding to expired US Patent No. 5,079,352. The purchase of this product includes a limited, non-transferable immunity from suit under the foregoing patent claims for using only this amount of product for the purchaser's own internal research. No right under any other patent claim and no right to perform commercial services of any kind, including without limitation reporting the results of purchaser's activities for a fee or other commercial consideration, is conveyed expressly, by implication, or by estoppel. This product is for research use only. Diagnostic uses under Roche patents require a separate license from Roche. Further information on purchasing licenses may be obtained by contacting the Director of Licensing, Applied Biosystems, 850 Lincoln Centre Drive, Foster City, California 94404, USA.

### Diamond Taq<sup>®</sup> is a registered trademark of Eurogentec S.A.

## For further information please contact our Customer Help Desk:

### For Europe:

e-mail: info@eurogentec.com

Tel.: +32 4 372 76 65

### For USA:

e-mail: info.anaspec@eurogentec.com

Tel.: +1 510 791 9560