

WHICH TUBES FOR YOUR BIORUPTOR® PLUS?



Catalog number	C30010013	C30010010	C30010012	C30010009	C20000010	14-432-22 (Fisher Scientific)
Recommended sample volume	50 µl - 100 µl	100 µl - 300 µl	500 µl - 2 ml	500 µl - 2 ml	500 µl - 2 ml	5 ml - 15 ml

APPLICATIONS	✓	✓	✓	✓	✓	✓
DNA shearing	✓					
Chromatin shearing	✓	✓	✓	✓		
FFPE DNA Extraction		✓				
Cell lysis	✓	✓	✓	✓	✓	✓
Bacterial cell lysis	✓	✓	✓	✓		
Tissue disruption for protein extraction		✓*			✓*	
Tissue disruption for RNA extraction						✓

✓ = validated tubes

* Protein extraction beads (C20000021) must be used



If you have any question, please contact the specialist at <https://www.diagenode.com/en/pages/support>